DAILY WEAR
(FREQUENT REPLACEMENT)
AND EXTENDED WEAR

PATIENT INSTRUCTION GUIDE

ACUVUE®
BRAND CONTACT LENSES

CE
2797

This patient instruction guide refers to the following ACUVUE® Brand Contact Lenses which are individually listed in **Table 1** below and are referenced as such unless stated otherwise.

Table 1	Intended use and wear schedule			
Lens type and Brand name	Daily wear (Frequent replacement)	Extended wear	Material	Packaging solution
ACUVUE® Brand Spherical Contact Lenses – Visibility Tinted with UV-Blocker				
ACUVUE® 2 Brand Contact Lenses	●	●	etafilcon A	①
ACUVUE OASYS® Brand Contact Lenses with HYDRACLEAR® PLUS	③ ●	●	senofilcon A	②
ACUVUE® VITA® Brand Contact Lenses	③ ●		senofilcon C	②
ACUVUE® OASYS Contact Lenses with Transitions™	③ ④ ●		senofilcon A	②
ACUVUE® Brand Contact Lenses for ASTIGMATISM – Visibility Tinted with UV-Blocker				
ACUVUE OASYS® Brand Contact Lenses for ASTIGMATISM with HYDRACLEAR® PLUS	③ ●	●	senofilcon A	②
ACUVUE® Brand Contact Lenses for PRESBYOPIA – Visibility Tinted with UV-Blocker				
ACUVUE OASYS® Brand Contact Lenses for PRESBYOPIA with HYDRACLEAR® PLUS	③ ●	●	senofilcon A	②

Packaging solution

- ① Borate buffered saline. ② Borate buffered saline with methyl ether cellulose.

Material content

- ③ Lens material contains silicone and meets Class 1 UV-absorbing standards with transmissibility of less than 1% UVB (280-315nm) and 10% UVA (316-380nm) radiation. All other ACUVUE® products meet Class 2 UV-absorbing standards with transmissibility of less than 5% UVB and 50% UVA radiation.
- ④ Lens material contains a photochromic¹ additive which dynamically absorbs visible light in the range from 380 nm to 780 nm to a minimum 84% transmittance in the inactivated (closed) state and to a minimum of 23% transmittance in the activated (open) state dependent on the lens thickness and the level of absorbed UV and high energy visible (HEV) radiation.

The following symbols may appear on the labels or cartons of your ACUVUE® Brand Contact Lenses.

Symbol	Definition	Symbol	Definition
	Caution, Consult instructions for use	MID / M	"Medium" near ADD
	Manufacturer	HGH / H	"High" near ADD
	Date of manufacture		Quality System certification symbol
	Use by date (Expiration Date)		UV-blocking
	Batch code		UV-blocking
	Sterilized using steam heat		Fee paid for waste management
DIA	Diameter		CAUTION: US federal law restricts this device to sale by or on the order of a licensed practitioner
BC	Base curve		Lens orientation correct
D	Dioptre (lens power)		Lens orientation incorrect (lens inside out)
CYL	Cylinder power		"Identification mark" for paper containers and wrapping
AXIS	Axis		"Identification mark" for composite materials
1-MONTH	1-Month replacement		Store away from direct sunlight
	Authorised Representative in the European Community		
	Do not use if package is damaged		
MAX ADD	Highest near addition that can be corrected		
LOW / L	"Low" near ADD		

¹Photochromic is defined as a reversible change in color or shade when exposed to specific types of light of sufficient intensity. In the absence of activating light, the color returns to clear.

Contents

Introduction.....	5
Intended use	6
Wear schedule	7
Contra-indications	8
Warnings – what you should know about contact lens wear	9
Adverse reactions	10
Recognising problems and what to do	11
A clean routine	12
The lens packaging.....	13
Inserting your lenses	14
Placing the lens on the eye	15
Centring the lens.....	16
Care for a sticking (non-moving) lens	17
Care for a dried out (dehydrated) lens.....	17
Removing your lenses	17
Caring for your lenses	18
Looking after your lens case.....	20
In the event of running out of lenses.....	20
Instructions for the presbyopic patient (monovision and multifocal)	21
Other important information / Precautions.....	22
Summary overview.....	24
Reporting adverse reactions	25
Your prescription.....	25

Introduction

As with any contact lens, ACUVUE® Brand Contact Lenses are medical devices. The ongoing support of your Eye Care Professional for regular eye checks is essential to the long-term protection of your health and vision.

For your eye health, it is important that your contact lenses be worn only as prescribed by your Eye Care Professional. He or she will review with you **all risks associated with contact lens wearing** and will provide instructions for lens handling and care, including how to safely and easily open the packaging. You will also be taught how to properly apply and remove lenses. This booklet will reinforce those instructions.

If you have any questions, always ask your Eye Care Professional.

Intended use

ACUVUE® Brand Spherical Contact Lenses are intended for Daily Wear and for the optical correction, myopia (short-sightedness) and hyperopia (long-sightedness) in persons with healthy eyes who have 1.00D or less of astigmatism.

ACUVUE® OASYS with Transitions™ are photochromic contact lenses also indicated for the attenuation of bright light as they contain a photochromic additive which dynamically absorbs visible light.

ACUVUE®2 and ACUVUE OASYS® Brand Contact Lenses with HYDRACLEAR® PLUS are also intended for Extended Wear.

ACUVUE® Brand Contact Lenses for ASTIGMATISM are intended for Daily Wear and for the optical correction of myopia (short-sightedness) and hyperopia (long-sightedness) in persons with healthy eyes who may have astigmatism.

ACUVUE OASYS® Brand Contact Lenses for ASTIGMATISM with HYDRACLEAR® PLUS are also intended for Daily Wear or Extended Wear.

ACUVUE® Brand Contact Lenses for PRESBYOPIA are intended for Daily Wear or Extended Wear for the optical correction of myopia (short-sightedness) and hyperopia (long-sightedness) in presbyopic persons with healthy eyes who have 0.75D or less of astigmatism.

ACUVUE OASYS® Brand Contact Lenses with HYDRACLEAR® PLUS are also indicated for therapeutic use as a bandage lens for certain ocular conditions.

Please consult your Eye Care Professional for further information. You should never self-treat with a contact lens or eye medications any condition without first being seen by your Eye Care Professional.

WARNING: All ACUVUE® Brand Contact Lenses have UV blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear such as UV-absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. You should continue to use UV-absorbing eyewear as directed by your Eye Care Professional.

Note:

Long term exposure to UV radiation is one of the risk factors associated with cataracts. Exposure is based on a number of factors such as environmental conditions (altitude, geography, cloud cover) and personal factors (extent and nature of outdoor activities). UV-blocking contact lenses help provide protection against harmful UV radiation. However, clinical studies have not been done to demonstrate that wearing UV-blocking contact lenses reduces the risk of developing cataracts or other eye disorders. Consult your Eye Care Professional for more information.

Wear schedule

Your Eye Care Professional should determine the appropriate wear and replacement schedule based upon your history and ocular examination.

Daily Wear – Frequent Replacement

ACUVUE® Brand Contact Lenses prescribed for Daily Wear – Frequent Replacement (less than 24 hours, while awake), are to be discarded and replaced every two weeks.

ACUVUE® VITA® Brand Contact Lenses prescribed for Daily Wear – Frequent Replacement (less than 24 hours, while awake), are to be discarded and replaced every month.

All ACUVUE® Brand Contact Lenses when prescribed for Daily Wear – Frequent Replacement are to be cleaned, rinsed and disinfected each time the lens is removed using a chemical disinfection system only.

Extended Wear

ACUVUE® Brand Contact Lenses prescribed for Extended Wear, (greater than 24 hours, including while asleep) as shown in **Table 1** may be used continuously for up to 7 days/6 nights and should be discarded upon removal. Your Eye Care Professional will determine if you are suitable for extended wear before lenses are worn in this way. When used in this way, no cleaning or disinfection is required.

It is recommended that you first be started on a daily wear schedule. If successful, then a gradual introduction of extended wear can be followed as determined by your Eye Care Professional.

Once removed, it is recommended that the lens remains out of the eye for a period of rest overnight or longer.

Contra-indications

When worn for vision correction, do not use ACUVUE® Brand Contact Lenses when you have any of the following conditions:

- Inflammation or infection in or around the eye or eyelids.
- Any eye disease, injury or abnormality that affects the corneas, conjunctiva or eyelids.
- Any previously diagnosed condition that makes contact lens wear uncomfortable.
- Severe dry eye.
- Reduced corneal sensitivity (corneal hypoesthesia).
- Any systemic disease that may affect the eye or be exaggerated by wearing contact lenses.
- Allergic reactions of ocular surfaces or surrounding tissues that may be induced or exaggerated by wearing or use of contact lens solutions.
- Allergy to any ingredient, such as mercury or Thimerosal, in a solution which is to be used to care for your contact lenses.

- Any active corneal infection (bacterial, fungal, protozoal or viral).
- If eyes become red or irritated.

For Therapeutic Use, your Eye Care Professional may prescribe ACUVUE OASYS® Brand Contact Lenses with HYDRACLEAR® PLUS to aid in the healing process of certain ocular conditions that may include those listed above.

Warnings – what you should know about contact lens wear

It is essential that you follow your Eye Care Professional's directions for the proper use of your contact lenses. You are advised of the following warnings pertaining to contact lens wear:

- Problems with contact lenses or lens care products could result in serious injury to the eye. You should be cautioned that proper use and care of contact lenses and lens care products, including lens cases, are essential for the safe use of these products.
- Eye problems, including corneal ulcers, can develop rapidly and lead to loss of vision.
- Studies have shown that the risk of eye problems, including ulcerative keratitis is greater for extended wear contact lens users than for daily wear users.
- When daily wear users wear their lenses overnight (outside the approved indication), the risk of eye problems, including ulcerative keratitis is greater than among those who do not wear them overnight†.
- The overall risk of eye problems, including ulcerative keratitis may be reduced by carefully following directions for lens care, including cleaning the lens case.
- Studies have shown that the risk of eye problems, including ulcerative keratitis among contact lens users who smoke is greater than among non-smokers.
- If you experience eye discomfort, excessive tearing, vision changes, redness of the eye or other problems, you should immediately remove your lenses and promptly contact your Eye Care Professional.

† New England Journal of Medicine, September 21, 1989, 321 (12), pp 773-783.

- It is recommended that you see your Eye Care Professional routinely as directed.
- **DO NOT** expose contact lenses to **water** during swimming, other water sports or bathing as this could increase the risk of serious eye infection from microorganisms which **could lead to vision loss**. If lenses have been submersed in water, the patient should discard and replace them with a new pair. The Eye Care Professional should be consulted for recommendations regarding wearing lenses during any activity involving water.
- Do not reuse or “top-off” old solution left in your lens case since solution reuse reduces effective lens disinfection and could lead to severe infection, vision loss, or blindness. “Topping-Off” is the addition of fresh solution to solution that has been sitting in your case.
- Discard any remaining solution after the recommended time period indicated on the bottle of multi-purpose solution used for disinfecting and soaking your contact lenses. The discard date refers to the time you can safely use the contact lens care product after the bottle has been opened. It is not the same as the expiration date, which is the last date that the product is still effective before it is opened.
- Using multi-purpose solution beyond the discard date could result in contamination of the solution and can lead to severe infection, vision loss, or blindness. To avoid contamination, DO NOT touch the tip of the container to any surface. Replace cap after using. DO NOT transfer to other bottles or containers.
- Do not store your lenses or rinse your lens case with water or any non-sterile solution. Only use fresh multi-purpose solution so you do not contaminate your lenses or lens case. Use of non-sterile solution can lead to severe infection, vision loss, or blindness.

Adverse reactions

Be aware that the following problems may occur when wearing contact lenses:

- Your eyes may burn, sting and/or itch.
- There may be less comfort than when the lens was first placed on the eye.
- There may be a feeling of something in your eye.

- There may be the potential for some temporary harm due to peripheral infiltrates, peripheral corneal ulcers and corneal erosion. There may be the potential for other physiological observations, such as local or generalised oedema, corneal neovascularisation, corneal staining, injection, tarsal abnormalities, iritis and conjunctivitis, some of which are clinically acceptable in low amounts.
- There may be excessive watering, unusual eye secretions or redness of your eye.
- Poor vision, blurred vision, rainbows or halos around objects, sensitivity to light (photophobia) or symptoms of eye dryness may also occur if your lenses are worn continuously or for too long a time.

During therapeutic use of ACUVUE OASYS® Brand Contact Lenses with HYDRACLEAR® PLUS, an adverse effect may be due to the original disease or injury may be due to the effects of wearing a contact lens. There is a possibility that the existing disease or condition might become worse when a soft contact lens for therapeutic use is used to treat an already diseased injured eye. To avoid serious eye damage, you should contact your Eye Care Professional IMMEDIATELY if there is an increase in symptoms while wearing the lens.

If any of the above symptoms occur, a serious eye condition such as infection, corneal ulcer, neovascularisation or iritis may be present. You should immediately be seen by your Eye Care Professional, so that the problem can be identified and treated, if necessary, in order to avoid serious eye damage.

Recognising problems and what to do

You should conduct a simple self-examination at least once a day.

Ask yourself:

1. How do the lenses feel on my eyes?
2. How do my eyes look?
3. Do I continue to see well?

If you notice any problems, you should IMMEDIATELY REMOVE YOUR LENS. If the problem or discomfort stops and the lens appears undamaged, clean and rinse the lens with a recommended soft contact lens care solution, and reinsert the lens. If after reinserting the lens, the problem continues, discard the lens and place a new fresh lens on the eye.

If after inserting the new lens, the problem continues, IMMEDIATELY REMOVE THE LENS AND CONTACT YOUR EYE CARE PROFESSIONAL.

REMEMBER – SYMPTOMS ARE WARNING SIGNS. IF IN DOUBT, TAKE THEM OUT.

A clean routine

Preparing the lens for insertion

Cleanliness is the first and most important aspect of proper contact lens care.

Establish a routine of good hygiene for handling your lenses.

- **Always** wash your hands thoroughly with warm water, a mild soap, rinse carefully and dry with a clean lint-free towel before touching your lenses to reduce the chance of getting an infection.
- Do not use cosmetics, soaps containing cold cream, lotions or creams before handling your lenses. It is best to insert your lenses before putting on make-up.
- Keep your eyes closed when using hairspray or other aerosols.
- Always follow the instructions in this booklet and any advice given to you by your Eye Care Professional for the correct handling, insertion, removal and use of your lenses.
- Never wear lenses longer than the period prescribed.

The lens packaging

To open the box, locate the opening flap on the front and pull up to break the seal.

Each lens is in its own package, designed specifically to maintain sterility while sealed. To close the box for storage, just tuck in the flap.

DO NOT USE if the sterile blister package is opened, damaged or after the expiry date shown.⌘

Always confirm the lens parameters (e.g. diameter (DIA), base curve (BC), lens power (D), etc.) printed on the multi-pack and on the individual lens package match your prescription. DO NOT use if there is a mismatch.

To open an individual lens package, follow these simple steps:

1. First separate one lens from the strip of lenses. Be careful not to break the seal on any of the other lens packages.
2. Shake the package to allow the lens to float free in the solution.
3. Peel back the foil. Occasionally, a lens may stick to the inside of the foil or to the package itself. This will not affect the sterility of the lens, which is still perfectly safe to use.
4. Handle your lenses carefully with your fingertips and be careful to avoid damage with your fingernails. It is helpful to keep these smooth and short.
5. Carefully remove the lens by sliding it up the side of the container. Never use tweezers or other tools.

Inserting your lenses

To avoid mix-ups, develop the habit of always inserting the first lens in your right eye.

Before inserting the lens, check to see that it is a single, moist, clean lens that is free of any nicks or tears. If it appears damaged, throw it away and use the next lens.

Check that the lens has not turned inside out.

Place the lens on the tip of your forefinger and check its profile. The lens should assume a natural, curved, bowl-like shape. If the lens edges tend to point outward, the lens is inside out. Another method is to gently squeeze the lens between the thumb and forefinger. The edges should turn inward. If the lens is inside out, the edges will turn slightly outward.

✓ Correct

✗ Incorrect

Or

Place the lens on the tip of your index finger and, looking up at the lens, locate the numbers 123. 1-2-3 indicates correct orientation while a reverse of 1-2-3 indicates the lens is inside out. If the lens is inside out (reverse 1-2-3), invert the lens and locate the numbers again to confirm correct lens orientation.

Note that the 1-2-3 mark is not present on ACUVUE® Brand Contact Lenses for ASTIGMATISM and ACUVUE OASYS® Brand Contact Lenses for PRESBYOPIA with HYDRACLEAR® PLUS.

Note:

If you need to rinse the lens before you insert it into your eye, use only fresh, sterile rinsing solution, as recommended by your Eye Care Professional.

NEVER USE TAP WATER.

Placing the lens on the eye

1. Remember to start with your right eye. Once the lens has been examined and you are sure it is not inside out, place it on the tip of your forefinger.
2. Place the middle finger of the same hand close to your lower eyelashes and pull down the lower lid.
3. Use the forefinger or middle finger of the other hand to lift the upper lid and place the lens on the eye.
4. Gently release both lids and blink.
5. Repeat these steps for the left lens.

There are other methods of lens placement. If the above method is difficult for you, your Eye Care Professional can provide an alternative.

Centring the lens

Usually, the lens centres itself automatically on the middle of your eye when you insert it and will very rarely be displaced onto the white of your eye during wear. However, this can occur if insertion and removal are not performed properly. To centre a lens, follow either of these methods:

- a. Close your eyelids and gently massage the lens into place through the closed lid.

Or

- b. Gently manipulate the off-centred lens onto the middle of your eye while the eye is opened, using finger pressure on the edge of the upper or lower lid.

Note:

If your vision is blurred after inserting the lens, check for the following:

- The lens may not be centred on the eye. Check the instructions above.
- If the lens is in the correct position, remove it and look for the following:
 - a. Cosmetics or oils on the lens. Clean and rinse with a recommended soft contact lens care solution, and if the cosmetic or oil has been removed, reinsert the lens. If your vision is still blurred, dispose of the lens and insert a new fresh lens.
 - b. The lens may be on the wrong eye.
 - c. The lens may be inside out, which would also make it less comfortable than normal.

If you find that your vision is still blurred after checking the above possibilities, remove the lens and consult your Eye Care Professional.

Care for a sticking (non-moving) lens

If a lens sticks on your eye, you may be prescribed a lubricating or wetting solution by your Eye Care Professional. In this way you can wet your lenses while you are wearing them, to make them more comfortable.

Just apply a few drops and wait until the lens begins to move freely on the eye. If this does not solve the problem, consult your Eye Care Professional **IMMEDIATELY.**

Care for a dried out (dehydrated) lens

If any ACUVUE® Brand Contact Lens is off the eye and exposed to air for a prolonged period, its surface may become dry and gradually become brittle. If this should occur, discard the lens and use a new one.

Removing your lenses

Always remove the same lens first.

Wash with soap & warm water, rinse and dry your hands thoroughly and follow the same rules for hygiene as already described in 'a clean routine'.

Caution:

Always ensure that the lens is on the middle of your eye before attempting to remove it.

You can determine this by covering the other eye and if vision is blurred, the lens is either on the white of the eye or it is not on the eye at all. To locate the lens, inspect the upper area of the eye by looking down into a mirror while pulling the upper lid up. Then inspect the lower area by pulling the lower lid down.

Once you have found the lens, you can remove it by using the Pinch Method or any other method recommended by your Eye Care Professional.

The Pinch Method

1. Look up and slide the lens down to the white of your eye using your forefinger.
2. Gently pinch the lens between your thumb and forefinger and remove the lens.

Caring for your lenses

For the continued safe and comfortable use of your lenses, it is important to follow the instructions given to you by your Eye Care Professional. If you remove your lenses and plan to wear them again, carefully read the instructions provided for proper cleaning, rinsing, disinfecting and storage.

Failure to follow the correct lens care regime may result in the development of serious eye problems, as described in the section entitled '**Warnings – what you should know about contact lens wear**'.

Cleaning and rinsing are essential to remove mucus, secretions and deposits which may have accumulated during use. Do this immediately after removing your lenses and prior to disinfection.

Harmful germs can only be removed by cleaning, rinsing and disinfecting.

To adequately disinfect the lenses, rub and rinse the lenses according to the recommended lens rubbing and rinsing times in the labelling of the multi-purpose solution.

Follow the instructions for use of the lens care system (chemical not heat) recommended by your Eye Care Professional.

Different products cannot always be used together and not all products are safe for use with all lenses. Do not alternate or mix lens care systems unless indicated on the lens care labelling.

When using hydrogen peroxide lens care systems, use **ONLY** the lens case provided with the hydrogen peroxide care system. This case is specially designed to neutralise the solution. Failure to use the specialised case will result in severe stinging, burning, and injury to the eye.

- **Always** wash with soap & warm water, rinse and dry hands before handling lenses.
- Use fresh, unexpired lens care solutions.
- **NEVER** change your solutions without consulting your Eye Care Professional.
- Thermal systems that use heat can damage your lenses.
- **NEVER** store your lenses in saline only. Saline will not protect the lenses from contamination and may result in an eye infection.
- **NEVER** use solutions recommended for rigid gas permeable (RPG) contact lenses.
- Always store individual unopened ACUVUE® OASYS with Transitions™ blisters out of direct sunlight.
- Always store worn ACUVUE® OASYS with Transitions™ contact lenses in the lens case and out of direct sunlight.
- Since some lens materials contain silicone, as shown in Table 1, the wettability may differ when different lens care products are used. Your Eye Care Professional should recommend a care system that is appropriate for ACUVUE® Brand Contact lenses.
- When you use sterile solutions that do not contain preservatives, they should be thrown out after the time specified in the directions.
- **NEVER** put lenses in your mouth or use anything other than the recommended solutions for lubricating or wetting your lenses.
- **NEVER** rinse them in tap water, since this can contain many impurities that can contaminate or damage your lenses and may lead to eye infection or injury.
- Clean the right lens first, to avoid mix-ups.

- Put each lens into the correct chamber of the lens storage system and make sure they are completely immersed in the storage solution when they are not being worn. If lenses are left out for long periods, they may dry out and become brittle. If this happens, throw them away.
- **NEVER** re-use the solution in your lens case.
- Seek the advice of your Eye Care Professional if your lenses are to be stored for extended periods.
- Lenses prescribed for reusable frequent replacement program should be thrown away after the recommended wearing period prescribed by your Eye Care Professional.

Looking after your lens case

Since lens cases can be a source of bacteria, after use they should be emptied, cleaned and rinsed with recommended sterile solutions and allowed to air dry. Your lens case should also be replaced regularly as advised by the lens case manufacturer or your Eye Care Professional.

NEVER USE TAP WATER TO RINSE YOUR LENS CASE.

In the event of running out of lenses

It is important to make sure you always have an adequate supply of replacement lenses.

To avoid running out, you will need to allow time to order and collect or to be supplied with your replacement lenses from your Eye Care Professional.

If you run out of lenses for any reason, you should wear your spectacles.

Although ACUVUE® Brand Contact Lenses maybe your preferred vision correction, spectacles are an essential back-up for all contact lens wearers.

Instructions for the presbyopic patient (monovision and multifocal)

If your Eye Care Professional explains that you require monovision correction, you should take note of the following:

As with any type of lens correction, there may be a visual compromise.

In some cases, monovision lenses can reduce visual acuity and depth perception for distance and near tasks. Some patients have difficulty adapting to this.

Symptoms such as mild blurring and variable vision may last briefly or for several weeks, as your eyes learn to adapt. The longer these symptoms persist, the poorer your chances for successful adaptation.

During this period, it is advisable to only wear these lenses in familiar situations that are not visually demanding. For example, until your eyes have adjusted, it would be wise to avoid driving.

Some patients may also need to wear spectacles over their lenses to provide the clearest vision for critical tasks.

- In some cases, patients will never be fully comfortable functioning in poor light when driving at night, for example. If this happens you may be prescribed additional lenses so that both eyes can be corrected when sharper distance vision is required.
- The decision to be fitted with monovision lenses should always be made after careful consultation with an Eye Care Professional.
- You should follow the advice you are given to help you adapt to monovision lenses and you should always discuss fully any concerns or problems you may have during and after the adaptation period.

Other important information

Precautions

- Before leaving your Eye Care Professional, ensure that you are able to put the lenses on and remove them or have someone else available who can do this for you.
- **DO NOT** touch your contact lenses with your fingers or hands if they are not completely clean because tiny lens scratches may occur, causing unclear vision and/or injury to your eye.
- Always handle lenses carefully and avoid dropping them.
- Remove your lenses immediately if your eyes become red or irritated.
- Inform all your doctor that you are a contact lens wearer.
- Always seek advice from your Eye Care Professional before using any eye drops.
- Certain medications, such as antihistamines, decongestants, diuretics, muscle relaxants, tranquillisers and treatments for travel sickness may cause dryness of the eye, increased lens awareness or blurred vision. If you experience a problem, you should seek proper medical advice.
- If chemicals are splashed into your eyes: **FLUSH EYES IMMEDIATELY WITH TAP WATER. CONTACT YOUR EYE CARE PROFESSIONAL WITHOUT DELAY OR GO TO THE CASUALTY DEPARTMENT OF YOUR NEAREST HOSPITAL.**
- Those who use oral contraceptives could develop changes in vision or lens tolerance. Your Eye Care Professional will advise you about this problem.
- Do not change your lens type (e.g. brand name, etc.) or parameters (e.g. diameter, base curve, lens power, etc.) without consulting your Eye Care Professional.
- As with any contact lens, regular eye check-ups are essential to maintain healthy vision.
- Never allow anyone else to wear your lenses. They have been prescribed to fit your eye and to correct your vision to the degree necessary. Sharing lenses greatly increases the chances of eye infections.

- For healthy comfortable vision, it is important that ACUVUE® Brand Contact Lenses are worn only as prescribed by your Eye Care Professional. This booklet will act as a reminder of these instructions.
- Your Eye Care Professional should be kept fully informed about your medical history and will recommend a lens and care system specific to your needs.
- When the replacement period prescribed by your Eye Care Professional is over, ACUVUE® Brand Contact Lenses should be discarded and replaced with a new sterile pair.
- Avoid all harmful or irritating vapours and fumes whilst wearing lenses.
- Always inform your employer of being a contact lens wearer. Some jobs may require the use of eye protection equipment or may require that you do not wear contact lenses.
- **Be aware** that wearing a darkened ACUVUE® OASYS with Transitions™ lens on only one eye is not recommended because it may cause disturbances in your ability to accurately judge depth and the motion of objects. It may also create a cosmetic concern.
- **Never** wear ACUVUE® OASYS with Transitions™ contact lenses as protection against artificial light sources, such as sun lamps, lasers, etc or stare directly at the sun or at an eclipse.
- **Never** stare directly at the sun or at an eclipse with or without ACUVUE® OASYS with Transitions™ contact lenses.
- **Be aware** that data is not available on the safety and performance of driving with ACUVUE® OASYS with Transitions™ contact lenses for individuals age 50 and older who may have crystalline lens opacities (i.e. cataract).
- If you have any questions, always ask your Eye Care Professional.

Summary overview

As with all types of contact lenses, there are basic guidelines which must be followed to protect and enhance your sight:

- Always follow the instructions given to you for safe lens wear. Regular check-ups are important to maintain peak performance and healthy vision.
- High standards of hygiene are essential for safe lens wear.
- **NEVER** wear lenses for longer than the prescribed period.
- **CLEANING, RINSING** and **DISINFECTING** are **ESSENTIAL**, each time the lens is removed to be re-inserted.
- **NEVER** re-use the solution in your lens case. Use fresh solution each time lenses are removed.
- If you experience a problem such as red/irritated eye or blurred vision, **IMMEDIATELY** remove your lenses and contact your Eye Care Professional.
- Always have a pair of spectacles available so that you are not tempted to wear your lenses when they should be removed. Whilst contact lenses may be your preferred correction, spectacles are an essential back-up.

As a wearer of ACUVUE® Brand Contact Lenses, you will soon recognise the advantages over spectacles or other contact lenses.

For maximum benefit and performance it is vital that you recognise how important it is to follow these guidelines.

You should only wear your lenses after following the advice given to you by your Eye Care Professional and that contained in this booklet.

Reporting adverse reactions

Any adverse reactions experienced whilst wearing ACUVUE® Brand Contact Lenses should be reported to your Eye Care Professional who will then report these direct to the manufacturer.

Your prescription

ACUVUE® Brand Contact Lenses.

Lens brand

Right lens:

Power and base curve

Left lens:

Power and base curve

YOUR LENSES SHOULD BE REPLACED AS RECOMMENDED BY YOUR EYE CARE PROFESSIONAL.

Care regime:

--

Should you have any queries or concerns call your Eye Care Professional.

Address and telephone number of your Eye Care Professional:

--

Regular visits to your Eye Care Professional are important for clear and healthy eyes.

Follow-up visits:

1	<table border="1"><tr><td></td></tr></table>		2	<table border="1"><tr><td></td></tr></table>		3	<table border="1"><tr><td></td></tr></table>	
4	<table border="1"><tr><td></td></tr></table>		5	<table border="1"><tr><td></td></tr></table>		6	<table border="1"><tr><td></td></tr></table>	

THIS PAGE IS MEANT TO BE KEPT BLANK.

THIS PAGE IS MEANT TO BE KEPT BLANK.

THIS PAGE IS MEANT TO BE KEPT BLANK.

Manufactured by:

Please refer to carton for site of manufacture.

USA

Johnson & Johnson Vision Care Inc.
7500 Centurion Parkway
Jacksonville, Florida, 32256
U.S.A.

IRELAND

Johnson & Johnson Vision Care (Ireland)
The National Technology Park
Limerick
Ireland

EU Authorised Representative:

AMO Ireland
Block B, Liffey Valley Office Campus
Quarryvale, Co.
Dublin
Ireland

www.acuvuearabia.com

PIG 000371

Daily Wear & Extended Wear

PIG (English)

PIG000371

©Johnson & Johnson (Middle East) Inc 2019.

ACUVUE®, ACUVUE OASYS®, HYDRACLEAR® and VITA®
are registered trademarks of Johnson & Johnson Medical Limited.

Transitions™ is a trademark of Transitions Optical, Inc., used under license by TOL and JJVCI.

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AMO Ireland
Block B, Liffey Valley Office Campus,
Quarryvale, Co.
Dublin,
Ireland

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©Johnson & Johnson (Middle East) Inc 2019.

ACUVUE®, ACUVUE OASYS®, HYDRACLEAR® and VITA®
are registered trademarks of Johnson & Johnson Medical Limited.

Transitions™ is a trademark of Transitions Optical, Inc., used under license by TOL and JVC.

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